

Elmasteam Accessories

Product Information and Prices



Contents

- Water canister	3
- Evaporating dish	4
- Water treatment	5
- Elma Steam Descaler	6
- Wall bracket	7
- Compressed air	8
- Medical Adapter Set	9
Set of different nozzle attachments	9
Nozzels and Tweezer	10
Dimension definition	11
- Practicals	13
Plastic Tweezer	13
Stainless steel Tweezer	14
Screen basket	15
Luer-Lock Adapter for Elmasteam 8 basic & Elmasteam 4 ^s	16
Nozzel	17

Product Information and Prices



Elmasteam Accessories - Water canister

- Integrated level control
- Automatic filling of ElmaSteam 8 steam jet units with integrated pump
- Uninterrupted filling - steam jet unit can remain pressurized



Order numbers and list prices


Order number	1070254
List price	237,00 €

The prices indicated are net prices ex works, excluding transport costs, insurance, customs duties and taxes.

Technical Data

Ext. Dimensions W / D / H (mm)	235 / 195 / 305
Capacity (l)	10,0
Hose length (m)	2,0

Dimension definition



Product contents

▪ Water canister including connection hose and level control
▪ Funnel

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasteam Accessories - Evaporating dish

- Condensate and soiling are collected during evaporation, keeping the workplace clean
- Easy to transport and flexible use in the cleaning process
- The individual parts are made of stainless steel and plastic, have a smooth surface and are easy to clean
- The user-friendly holder for the Elmasteam handpiece can be positioned for left or right-handed use
- Bayonet openings in the base plate for attaching e.g. Elma pins, base comb
- The base tray can hold up to 2 liters of liquid



Order numbers and list prices

Order number	1127588
List price	260,00 €

The prices indicated are net prices ex works, excluding transport costs, insurance, customs duties and taxes.

Technical Data

Ext. Dimension W / D / H (mm)	455 / 385 / 285
Int. Dimension W ₁ / D ₁ (mm)	410 / 335
Weight (kg)	3,4


Product contents

▪ Bottom tray
▪ Grand plate with handpiece holder
▪ Splash guard with edge protection
▪ Elma Pin

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Dimension definition



Elmasteam Accessories - Water treatment

Water softening unit Brita purity C150 Quell ST

- reduces carbonate hardness in water
- decarbonizes water in four stages: 1. pre-filtration - 2. carbonate hardness reduction - 3. activated carbon filtration - 4. fine filtration
- PURITY C 30 % filter head with fixed cut-off Capacity1 at a carbonate hardness of 10 °dH = 2,086 liters
- positioning of the filter cartridge: horizontal or vertical
- maximum operating pressure: 8.6 bar



Order numbers and list prices

	Water softening unit	Filter Cartridge	filter head	Connection set	Reduction with seal (2x)	Hose for water connection	Extension hose
Order number	1084664	1084666	1084667	1084680	1084668	1059056	1065691
List price	290,00 €	85,00 €	35,00 €	85,00 €	18,00 €	18,50 €	50,00 €

The prices indicated are net prices ex works, excluding transport costs, insurance, customs duties and taxes.

Technical Data

	Water softening unit	Filter Cartridge	filter head	Connection set	Reduction with seal (2x)	Hose for water connection	Extension hose
Dimension W / D / H (mm)	117 / 104 / 419	104 / 104 / 410	--				
Dimensions	--				3/4" - 3/8"	Ø 10 mm; L= 1 m; 3/4" - 3/4" 90°	Ø 10 mm; L = 2 m; IG 3/4" - AG 3/4" 90°
Capacity at 20° german hardness (L)	500		--				
							

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Product Information and Prices



Elmasteam Accessories - Elma Steam Descaler

- Descaler, citric acid-based, for Elmasteam 8
- 1 packet is easy and safe to dose
- Odorless and tasteless
- Removes and prevents the build-up of limescale and boiler scale
- Descaling extends the service life, saves energy and expensive repairs



Order numbers and list prices

Order number	8000878
List price	12,00 €

The prices indicated are net prices ex works, excluding transport costs, insurance, customs duties and taxes.

Product contents

- 5 packets of Elma Steam Descaler (45g each)
- Brief instructions descaling


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Product Information and Prices



Elmasteam Accessories - Wall bracket

- Parts are made of stainless steel and plastic
- Mounting set for fixing the Elmasteam 4^s to the wall




Order numbers and list prices

Order number	1058677
List price	35,00 €

The prices indicated are net prices ex works, excluding transport costs, insurance, customs duties and taxes.

Wall bracket



Product contents

▪ 2x Foot (1)
▪ 2x Spacer (2)
▪ 2x Dowel (3)
▪ 2x Screw (4)

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Product Information and Prices



Elmasteam Accessories - Compressed air

- the maximum inlet pressure of 6 bar can be reduced accordingly to below 6 bar with the pressure reducer on the Elmasteam 8 basic and the Elmasteam 8 med



Order numbers and list prices


	Pressure reducer	Adapter set compressed air
Order number	1092886	1080066
List price	130,00 €	19,50 €

The prices indicated are net prices ex works, excluding transport costs, insurance, customs duties and taxes.

Technical Data

	Pressure reducer	Adapter set compressed air
Ext. Dimension W / D / H (mm)	60 / 75 / 125	--

Dimension definition



Product contents

Druckminderer	Adapterset Druckluft
■ 1x Pressure reducer including 2x push-in fitting and 1x 1m hose	■ 1x Plug connector for hose outer diameter 6mm to 8mm
	■ 1x Plug connector for hose outer diameter 6mm to 10mm
	■ 2x 0,35m Hose with outer diameter 6mm

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasteam Accessories - Medical Adapter Set

Set of different nozzle attachments

- connectable to an Elmasteam 8 med steam cleaner with standardized Luer-Lock adapter
- allows the cleaning of components and instruments with difficult-to-access areas
- storage rack for 8 different attachments
- reliable cleaning for a wide range of applications

Order numbers and list prices

Order number	1028820
List price	600,00 €

The prices indicated are net prices ex works, excluding transport costs, insurance, customs duties and taxes.

Technical Data

Ext. Dimension W / D / H (mm)	130 / 70 / 560
Height rack H ₁ (mm)	40


Product contents

■ 1x Storage rack	■ 2x Catheter attachment
■ 2x MIS hollow needle	■ 1x Hose connection
■ 2x Hollow needle attachment	■ 1x Pliers
■ 1x Luer-Lock adapter	

The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Dimension definition



Elmasteam Accessories - Medical Adapter Set










Nozzels and Tweezer

Order numbers and list prices

	MIS hollow needle	MIS hollow needle with brush	Hollow need- le attachment straight	Hollow need- le attachment circumflex	Catheter attachment large	Catheter attachment medium	Hose connection large	Adapter Luer-Lock female / female	Pliers
Order number	1028822	1055030	1028823	1028824	1028825	1028826	1028827	1028828	1028831
List price	199,00 €	178,00 €	36,50 €	36,50 €	40,00 €	40,00 €	31,50 €	25,50 €	48,50 €

The prices indicated are net prices ex works, excluding transport costs, insurance, customs duties and taxes.

Technical Data

	MIS hollow needle	MIS hollow needle with brush	Hollow need- le attachment straight	Hollow need- le attachment circumflex	Catheter attachment large	Catheter attachment medium	Hose connection large	Adapter Luer-Lock female / female	Pliers
Total length L ₁ (mm)	525	338	103	88	44	44	26	21	115
Usable length L ₂ (mm)	501	315	84	61	29	29	12	9	16
Diameter D (mm)	2,5	2,5	2	2	2-9	2,5-6,5	3,5	9	--
Width W (mm)	--	--	--	--	--	--	--	--	29
Height H (mm)	--	--	--	--	--	--	--	--	15
									

The drawings are meant as a visualization aid and are not true to scale.


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasteam Accessories - Medical Adapter Set

Dimension definition


MIS hollow needle

for tube shafts,
lateral steam outlet




MIS hollow needle

with brush, for tube shafts,
steam outlet behind brush hem




Hollow needle attachment

straight, front steam outlet



Hollow needle attachment

curved, front steam outlet




The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.


Elmasteam Accessories - Medical Adapter Set

Dimension definition


Catheter attachment
large,
front steam outlet




Catheter attachment
medium,
front steam outlet




Hose connection
large,
front steam outlet



Adapter Luer-Lock
fem. / fem.,
front steam outlet



Pliers




The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasteam Accessories - Practicals

Plastic Tweezer

- for precise gripping of cleaning goods without contact of the hands with the cleaning agent



Order numbers and list prices



Order number	2000001536
List price	3,50 €

The prices indicated are net prices ex works, excluding transport costs, insurance, customs duties and taxes.

Technical Data

Ext. Dimension W / D / H (mm)	25 / 130 / 10

Dimension definition



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasteam Accessories - Practicals

Stainless steel Tweezer

- the front part of the tweezer has a coating of PVC
- for precise gripping of cleaning goods without contact of the hands with the cleaning agent
- Tweezer 1118673 have an integrated lock that secures the tweezer when closed






Order numbers and list prices


	Straight short	Straight long	Curved with lock
Order number	376071000	1118674	1118673
List price	14,00 €	22,00 €	16,50 €

The prices indicated are net prices ex works, excluding transport costs, insurance, customs duties and taxes.

Technical Data

	Straight short	Straight long	Curved with lock
Ext. Dimension W / D / H (mm)	200 / 25 / 15	300 / 30 / 15	150 / 20 / 15
			

Dimension definition



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

Elmasteam Accessories - Practicals

Screen basket

- Stainless steel - PVC-coated
- Hinged
- Autoclavable



Order numbers and list prices


	Small	Big
Order number	3760720000	3760730000
List price	9,30 €	14,00 €

The prices indicated are net prices ex works, excluding transport costs, insurance, customs duties and taxes.

Technical Data

	Small	Big
Ext. Dimension W / D / H (mm)	155 / 45 / 40	185 / 75 / 70
Internal diameter Dia (mm)	35	64
Mesh size Y x Z (mm)	1,0 x 0,2	

Dimension definition



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.

Elmasteam Accessories - Practicals

Luer-Lock Adapter for Elmasteam 8 basic & Elmasteam 4⁵

- For connecting accessories with Luer-Lock connection, e.g. Medical Adapter Set



Order numbers and list prices


Order number	1091332
List price	69,00 €

The prices indicated are net prices ex works, excluding transport costs, insurance, customs duties and taxes.

Technical Data

Total length L_1 (mm)	33
Usable length L_2 (mm)	30
Diameter D (mm)	10

Dimension definition



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasteam Accessories - Practicals

Nozzel

- made of stainless steel and plastic
- control via foot switch allows cleaning with both hands



Order numbers and list prices



	Fixed Nozzle	Nozzle for Elmasteam basic handpiece
Order number	1076042	1076043
List price	69,00 €	20,00 €

The prices indicated are net prices ex works, excluding transport costs, insurance, customs duties and taxes.

Technical Data

	Fixed Nozzle	Nozzle for Elmasteam basic handpiece
Ext. Dimension W / D / H (mm)	140 / 20 / 75	45 / 10
Internal diameter Dia (mm)	2	
		

Dimension definition



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.